SILVERCEL® NON-ADHERENT

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EN Product Description

SILVERCEL® Non-Adherent Hydro-Alginate Antimicrobial Dressing/Packing with Silver is a sterile, non-woven pad composed of a high G (guluronic acid) alginate, carboxymethyl cellulose (CMC) and silver coated fibres, laminated to a perforated, non-adherent ethylene methyl acrylate (EMA) wound contact layer. The unique composition of the dressing manages exudates in moderate to heavily exuding wounds, which creates a favourable moist wound healing environment for effective wound management, and allows intact removal. The silver fibres kill a broad spectrum of microorganisms associated with the bacterial colonization and infection of wounds. Odour reduction results from the antibacterial effect in the dressing.

Indications

SILVERCEL® Non-Adherent Hydro-Alginate Antimicrobial Dressing is intended for use in the management of all moderate to heavily exuding partial and full-thickness chronic and acute wounds including: decubitus (pressure) ulcers, venous ulcers, diabetic ulcers, donor sites, traumatic and surgical wounds. As the product contains alginate it may assist in supporting the control of minor bleeding in superficial wounds. It is also suitable for use, under medical supervision, in the management of infected wounds, or wounds in which there is an increased risk of infection.

Contra-indications

SILVERCEL® Non-Adherent Hydro-Alginate Antimicrobial Dressing is not indicated for use on the following:

- Patients with a known sensitivity to alginates, ethylene methyl acrylate (EMA) or silver.
- Surgical implantation.
- Pregnant or lactating women due to the absence of specific information.

Precautions

- SILVERCEL® Non-Adherent Hydro-Alginate Antimicrobial Dressing is not intended to control heavy bleeding.
- As wound conditions improve and exudate levels decrease, it may be preferable to switch to a more appropriate dressing or moisten the dressing with normal saline solution prior to application.
- The dressing must be removed prior to patients undergoing Magnetic Resonance Imaging (MRI) examinations.
- Avoid contact with electrodes or conductive gels during electronic measurements, e.g. electrocardiograms (ECG) and electroencephalograms (EEG).
- Clinicians/Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.

Systemic antimicrobial therapy should be considered when wound infection is evident. SILVERCEL® Non-Adherent Hydro-Alginate Antimicrobial Dressing may be used, under medical supervision, in conjunction with systemic antibiotics.

DIRECTIONS FOR USE

Site Preparation

- Debride when necessary and irrigate the wound site in accordance with standard protocols.
- Remove excess solution from surrounding skin.

Dressing Selection

Select a size of SILVERCEL® Non-Adherent Hydro-Alginate Antimicrobial Dressing that is slightly larger than the wound.

Dressing Application

1. Cut (using clean scissors) or fold the dressing to fit the wound. Loosely pack deep wounds, ensuring the dressing does not overlap the wound margins. As wound conditions improve and exudate levels decrease, moisten the dressing with normal saline solution prior to application.
2. Cover and secure SILVERCEL® Non-Adherent Hydro-Alginate Antimicrobial Dressing with a non-occlusive secondary dressing such as one of the TIELLE® Hydropolymer Dressing range, as appropriate.

Dressing Change and Removal

1. Dressing change frequency will depend on wound condition and the level of exudate. Initially it may be necessary to change the dressing every 24 hours.
2. Reapply SILVERCEL® Non-Adherent Hydro-Alginate Antimicrobial Dressing when the secondary dressing has reached its absorbent capacity or whenever good wound care practice dictates that the dressing should be changed.
3. Gently remove the secondary dressing.
4. Gently remove the dressing from the wound bed and discard. If the wound appears dry, saturate the dressing with sterile saline solution prior to removal.
5. Irrigate the wound site prior to application of a new dressing.

Do not use if package is damaged.

The use by date of this product is printed on the packaging.

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©2018 Systagenix Wound Management Limited,
X-STATIC® is a registered trademark of -
Noble Fiber Technologies, Inc.,
Advanced Medical Solutions Limited,