




# ACTISORB™

## SILVER 220

ANTIMICROBIAL BINDING DRESSING





## **EN** Product Description

ACTISORB™ Silver 220 Antimicrobial Binding Dressing is a sterile primary dressing, comprised of activated charcoal cloth, impregnated with silver, within a spun bonded nylon envelope. Within the dressing, there is 220mg silver per 100g activated charcoal cloth equating to 33µg of silver per square centimeter of cloth.

### **Indications**

ACTISORB™ Silver 220 Dressing provides an effective barrier to bacterial penetration and for adsorbing offending odor resulting from wounds; the binding properties of the dressing trap bacteria, bacterial toxins and odor. ACTISORB™ Silver 220 Dressing may help reduce infection in partial and full-thickness wounds, including:

- pressure ulcers
- venous ulcers
- diabetic ulcers
- first and second-degree burns
- donor sites
- surgical wounds

ACTISORB™ Silver 220 Dressing is suitable for use under compression bandaging.

### **Contraindications**

ACTISORB™ Silver 220 Dressing should not be used on third-degree burns.

### **Warnings**

- Should not be used on patients with a known sensitivity to silver.
- For external use only.

### **Precautions**

- This product should not be cut, otherwise particles of activated charcoal may get into the wound and cause discoloration.
- The dressing performance may be impaired by excess use of petrolatum based ointments.
- Clinicians/Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.

## **DIRECTIONS FOR USE**

### **Site Preparation**

Prepare wound bed according to wound management protocol, ensuring removal of necrotic tissue.

## Dressing Application

- The properties of the dressing are most fully utilized when ACTISORB™ Silver 220 Antimicrobial Binding Dressing makes direct contact with the wound (see photo 2).
- Either side of ACTISORB™ Silver 220 Dressing can be applied to the wound bed.
- Do not cut dressing. Fold or pack as necessary.
- ACTISORB™ Silver 220 Dressing can be easily packed into deep wounds (see photo 3).
- The outer nylon envelope enables removal from the majority of wounds with minimal adherence or trauma.
- For lightly exuding wounds, ACTISORB™ Silver 220 Dressing can be used in combination with a non-adherent wound contact layer (see photo 4).
- Depending on the level of exudate, an absorbent secondary dressing may be placed on top of ACTISORB™ Silver 220 Dressing (see photo 5).
- Secure the dressing in a manner appropriate for the indication.

## Dressing Change and Removal

ACTISORB™ Silver 220 Antimicrobial Binding Dressing can remain in situ up to 7 days, dependent on the level of exudate, while the secondary absorbent dressing is changed as required. Initially it may be necessary to change ACTISORB™ Silver 220 Dressing every 24 hours (see photo 6) to ensure optimal dressing performance.

## STORAGE INFORMATION

Store at room temperature.

## HOW SUPPLIED

ACTISORB™ Silver 220 Dressing

- is supplied sterile
- package contains a single dressing
- should not be cut

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a properly licensed healthcare practitioner.

*This caution is not applicable outside the U.S.*

**Do not re-use (single use).**

**Do not resterilize.**

**Do not use if package is damaged.**

**The use by date of this product is printed on the packaging.**

**EN** SYMBOLS USED ON LABELLING



= Do not use if package is damaged.



= Do not re-use (single use).



= Do not re-sterilize.



= Sterilized using irradiation.



= Consult instructions for use.



= Caution.



= Batch number.



= Use by.



= Catalogue number.



= Manufacturer.



518044K.a



Manufacturer

©2018 Systagenix Wound Management Limited,  
Gargrave, North Yorkshire, BD23 3RX, UK



2018-07