 Let's Protect®

ACTISORB® PLUS 25

ACTIVATED CHARCOAL DRESSING WITH SILVER

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



 Systagenix

EN Product Description

ACTISORB® Plus 25 Activated Charcoal Dressing with Silver is a dressing composed of pure activated carbon impregnated with silver (25 µg/cm²). The dressing creates a favourable environment for effective wound healing by binding and immobilising micro-organisms which contaminate and infect wounds. The silver impregnation combats wound micro-organisms in the dressing, which reduces the bacterial colonisation of the wound and inhibits infection. ACTISORB® Plus 25 dressing is tolerated very well and because of its purely physical wound cleansing principle there is no fear of side effects or the development of resistance such as can occur when using antibiotics.

By eliminating offensive odours ACTISORB® Plus 25 dressing improves the quality of life for the patient and nursing personnel.

Indications

ACTISORB® Plus 25 dressing is suitable as the first therapeutic step in the management of all chronic wounds. It is indicated for fungating carcinomas, ulcerative, traumatic and surgical wounds where bacterial contamination, infection or odour occurs.

Precautions

Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.

DIRECTIONS FOR USE

Preparation of wound site

Dry, necrotic tissue must first be removed by surgical, enzymatic or autolytic debridement.

Application

- Prior to application, ACTISORB® Plus 25 dressing can be impregnated with saline solution or sterile water.
- The adsorptive properties of the dressing are most fully utilised when ACTISORB® Plus 25 dressing makes direct contact with the wound.
- Either side of ACTISORB® Plus 25 dressing can be placed next to the wound site.

- ACTISORB® Plus 25 dressing can be easily packed into deep wounds.
- The outer nylon surface enables removal from the majority of wounds without adherence or trauma.
- If required, ACTISORB® Plus 25 dressing can be placed on top of a non-adherent wound contact layer. In this case it is important to ensure that the mechanism of action of ACTISORB® Plus 25 dressing is not impaired by excess ointment or grease of the contact layer.
- Depending on the level of exudate an absorbent, secondary dressing may be placed on top of ACTISORB® Plus 25 dressing.
- Fix the dressing in place in a manner appropriate for the indication.
- Naturally, work according to a sterile procedure.

Changing the dressing

ACTISORB® Plus 25 dressing can remain in situ up to 7 days, dependent on the level of exudate, while the secondary absorbent dressing is changed as required. Initially it may be necessary to change ACTISORB® Plus 25 dressing every 24 hours.

N.B.: ACTISORB® Plus 25 dressing should not be cut, otherwise particles of activated carbon may get into the wound and cause discoloration.

Do not re-use.

Do not re-sterilize.

Do not use if package is damaged.

The use by date of this product is printed on the packaging.

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©2018 Systagenix Wound Management Limited,