
PROMOGRAN™

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 **KCI™**
AN ACELITY COMPANY

EN Product description

PROMOGRAN™ Matrix Wound Dressing is an advanced wound care device comprised of a sterile, freeze dried composite of 45% oxidized regenerated cellulose (ORC) and 55% collagen. In the presence of exudate the PROMOGRAN™ Dressing transforms into a soft, conformable, biodegradable gel, and thus allows contact with all areas of the wound. PROMOGRAN™ Dressing maintains a physiologically moist microenvironment at the wound surface. This environment is conducive to granulation tissue formation, epithelialization and rapid wound healing. PROMOGRAN™ Dressing is a primary dressing that can be cut to fit wound with sterile scissors and used in combination with either a semi-occlusive or non-occlusive secondary dressing. Prior to application in dry wounds saline solution should be used to hydrate PROMOGRAN™ Dressing.

Indications for use

PROMOGRAN™ Dressing is intended for the management of exuding wounds including:

- Diabetic ulcers
- Venous ulcers
- Pressure ulcers
- Ulcers caused by mixed vascular etiologies
- Full-thickness & partial thickness wounds
- Donor sites and other bleeding surface wounds
- Abrasions
- Traumatic wounds healing by secondary intention
- Dehisced surgical wounds

PROMOGRAN™ Dressing may be used under compression therapy with healthcare professional supervision.

Precautions

PROMOGRAN™ Dressing may be used when visible signs of infection are present in the wound area only when proper medical treatment addresses the underlying cause.

Contraindications

PROMOGRAN™ Dressing is not indicated for wounds with active vasculitis, third-degree burns, or patients with known sensitivity to ORC or collagen.

DIRECTIONS FOR USE

Wound Bed Preparation

Prepare wound bed per your standard wound care protocol and debride when necessary. Wound bed should be clear of any visual signs of infection.

Application

- For optimal effect, apply PROMOGRAN™ Dressing directly to the whole wound bed.
- For a wound with low or no exudate apply PROMOGRAN™ Dressing and hydrate with saline solution. This will initiate the transformation of the PROMOGRAN™ Dressing into a gel matrix.
- After hydration, through exposure to wound exudate or saline, the PROMOGRAN™ gel matrix will intimately come into contact with the wound surface.
- The biodegradable PROMOGRAN™ gel matrix is naturally absorbed into the body over time.
- In order to maintain a moist wound healing environment PROMOGRAN™ Dressing must be covered with a semi-occlusive dressing (e.g. BIOCLUSIVE™ PLUS Transparent Dressing, TIELLE® Hydropolymer Dressing or TIELLE® Plus Hydropolymer Dressing) or a non-occlusive secondary dressing and fixed to the skin with a nonirritating tape.
- After initial application, reapply PROMOGRAN™ Dressing to the wound up to every 72 hours depending upon the amount of exudate. It is not necessary to remove any residual PROMOGRAN™ Dressing during dressing changes.

Adverse Reactions

PROMOGRAN™ Dressing is contraindicated for patients with known hypersensitivity to the components of this product, i.e. ORC and Collagen. Discontinue use if signs of sensitivity appear.

Do not use if package is damaged.

Do not re-use.

Do not resterilize.

The use by date of this product is printed on the packaging.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a properly licensed healthcare practitioner.

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©2018 Systagenix Wound Management Limited,