# SOLUTION CONTAINER HANGER ARM EXTENSION

## INSTRUCTIONS FOR USE


Only for use with the V.A.C. ULTA™ Therapy Unit.

This component is intended to increase the height of the Solution Container Hanger Arm when using V.A.C. VERAFLU™ Therapy. For more information about the V.A.C. ULTA™ Therapy Unit, including detailed operating instructions, please consult the V.A.C. ULTA™ Therapy Unit User Manual.

## SETTING UP THE SOLUTION CONTAINER HANGER ARM EXTENSION ON THE V.A.C. ULTA™ THERAPY UNIT

## STORING THE SOLUTION CONTAINER HANGER ARM EXTENSION WHEN V.A.C. VERAFLU™ THERAPY IS NOT IN USE

## CARE AND CLEANING

This component may be reused following cleaning and disinfecting per hospital protocol. KCI recommends the Solution Container Hanger Arm Extension to be cleaned and disinfected:

- If it becomes soiled during patient use
- At least weekly
- Between patient use

## WARNINGS

Indications, Contraindications, Warnings, Precautions and other important Safety Information are contained in the V.A.C.ULTA™ Negative Pressure Wound Therapy System Safety Information included with the V.A.C.ULTA™ Therapy Unit.

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KCI USA, Inc.