 Let's Comfort™

---

# TIELLE®

HYDROPOLYMER ADHESIVE DRESSING  
WITH

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 Systagenix

## **EN** Product Description

TIELLE® Hydropolymer Adhesive Dressing provides an exudate handling system for low to moderately exuding wounds. It is a highly absorbent dressing. The absorbent material is a synthetic polymer. The island dressing maintains a moist environment that aids in the wound healing process and allows granulation to proceed under optimum conditions. During use the absorbent island gently expands as it takes up exudate. It allows the wound to remain moist thus encouraging autolytic debridement which may initially increase lesion size. This is normal and to be expected prior to wound granulation.

### **Indications**

TIELLE® Dressing is indicated for the management of low to moderately exuding wounds.

TIELLE® Dressing should be used under health care professional direction for the following indications:

- Pressure ulcers
- Lower extremity ulcers
  - Venous
  - Arterial
  - Mixed aetiology
- Diabetic ulcers
- Donor sites

TIELLE® Dressing is suitable for use under compression bandaging.

### **Precautions**

TIELLE® Dressing is not indicated for use on the following:

- Third-degree burns
- Lesions with active vasculitis as this type of ulcer needs more frequent observations by a healthcare professional

TIELLE® Dressings may be used when visible signs of infection are present in the wound area only when proper medical treatment

addresses the underlying cause.

## **DIRECTIONS FOR USE**

### **Site Preparation**

1. Prepare the wound according to wound management protocol.
2. Ensure skin surrounding the wound is dry.

### **Dressing Application**

**The size of the dressing selected should allow the absorbent island to overlap the wound edge by approximately 1 cm.**

1. Peel open the package and remove the dressing.
2. Partially peel back side backing papers. Position absorbent island centrally over wound site and smooth in place.
3. One at a time, peel away side backing papers while smoothing adhesive border onto intact skin.

### **Dressing Change and Removal**

It is essential for proper wound management to avoid trauma to the intact skin.

1. Change dressing when wound fluid is present at the edges of the foam pad. Do not allow exudate to accumulate under the backing. The dressing may be left in place up to 7 days depending upon the amount of exudate.
2. On removal lift one corner and carefully peel back. On fragile or friable skin, water or saline may be used to break the adhesive seal.

**Do not use if package is damaged.**

**Do not reuse (single use).**

**Do not resterilize.**

**The use by date of this product is printed on the packaging.**

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a properly licensed healthcare practitioner.

*This caution is not applicable outside the U.S.*

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©2018 Systagenix Wound Management Limited,